



AMPLATZER™ PFO OCCLUDER

**A LANDMARK DEVICE.
AND A TURNING POINT
FOR PFO CLOSURE.**



Abbott

A PIONEERING DEVICE WORLDWIDE

As the device that created the category, AMPLATZER™ PFO occluder has sustained leadership over decades of use by pursuing clinical evidence—even beyond an initial study end date—to become the first device supported by positive PFO trial results¹. Today, we continue to innovate around advancing patient safety and reducing risk for patients around the world.

A LANDMARK DEVICE. AND A TURNING POINT FOR PFO CLOSURE.

- Industry-leading device, designed for ease of use and effective closure
- Backed by the largest trial with the most extensive patient follow-up
- Demonstrated excellent safety and efficacy¹



WHY IS AMPLATZER PFO OCCLUDER RELIED UPON BY THOUSANDS OF PHYSICIANS AROUND THE WORLD?

GLOBALLY OVER
100K
IMPLANTS

WE SET THE STANDARD

- Pioneered treatment with a **PFO-specific** device
- **Over 100,000** devices implanted globally

5,810
PATIENT YEARS
OF DATA

WE RAISE THE BAR: THE LANDMARK RESPECT TRIAL¹

- Most extensive patient **follow-up**, **>2X more** than other PFO trials
- Only trial to include patients on anticoagulation therapy, a **real-world cross-section** of patients

990
PATIENTS IMPLANTED
WITH DEVICE
IN RCTS*

WE DEMONSTRATE EXCELLENCE

- **ZERO** device erosions, thrombus, or embolization events in **SIX**** published trials with **990** patients¹⁻⁶
- **94.2%** effective closure rates at 6 months¹

* RCTs=Randomized Clinical Trials

**Patients in device group of each trial implanted with AMPLATZER PFO occluder device: RESPECT = 465, PREMIUM = 119, PC = 191, CLOSE = 121, DEFENSE = 53, PRIMA = 41. RCTs=Randomized Clinical Trials

OFTEN IMITATED, NEVER MATCHED

Industry-leading device. Developed specifically for the treatment of PFO closure.

DESIGN ADVANTAGES THAT MAKE THE DIFFERENCE


DURABLE NITINOL WIRE MESH WITH POLYESTER FABRIC THREAD

Excellent visibility under fluoro

ASYMMETRIC DOUBLE DISC DESIGN

Minimizes material in the Left Atrium*

94%
CLOSURE RATE†
at 6 months in
RESPECT trial¹



INTAGLIO™ WIRE TREATMENT


Designed to reduce nickel leaching⁷

*On commonly used sizes (25mm and 35mm devices)

†Effective Closure

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

ONE STEP AHEAD: MINIMIZING COMPLEXITY IN CLOSURE



UNIQUE SELF-EXPANDING DISCS

Linked by a short-connecting waist, the discs align to the PFO without an additional “locking” step.

LOW PROFILE DELIVERY

8 F & 9 F introducer sheaths for multiple sizes (18, 25, 35 mm) enables treatment of patients with smaller vasculature.



FULLY RECAPTURABLE AND REPOSITIONABLE DESIGN

Allows confirmation of device placement prior to final release of the device.

THE TURNING POINT FOR PFO CLOSURE

Three trials published concurrently in the NEJM provide conclusive evidence of the superiority of PFO closure versus medical management in reducing risk of recurrent stroke.

	RESPECT ¹	REDUCE ⁸	CLOSE ⁴
Devices Used	100% AMPLATZER™ PFO Occluder	39% GORE [†] HELEX, 61% GORE Cardioform	51% AMPLATZER PFO Occluder; 49% other approved PFO devices
Patients	980	664	473
Follow-Up-Patient Years	5,810 (median 5.9 yrs)	2,232 (median 3.2 yrs)	NR* (mean 5.4 yrs)
Anticoagulant Allowed in Control Group?	Yes	No	No
	↓	↓	↓
Relative Risk Reduction	62% (Recurrent ischemic stroke of unknown mechanism)	77% (Recurrent ischemic stroke)	97% (Recurrent ischemic stroke)
Effective Closure	94.2% Freedom from >9 bubbles (Evaluated after 6 months)	94.5% Freedom from >25 bubbles (Evaluated after 12 months)	NR*

*Not reported

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

PUBLISHED DATA HIGHLIGHTS EXCELLENT SAFETY



LOW RISK OF DEVICE-RELATED EVENTS

	RESPECT ¹	REDUCE ⁸
Device Embolization/ Dislocation	0	2
Aortic Erosion/ Dissection	0	1
Device Thrombus	0	2

CLOSE Trial data not included due to device- and procedure-related events reported in combination.



LOW RISK OF ATRIAL FIBRILLATION (AF)

RATE (PER 100 PT YRS)	RESPECT ¹	REDUCE ^{8*}
Serious AF ^{**} / Flutter	0.22	0.65
Any AF/ Flutter	0.76	1.90

CLOSE Trial data not included as follow-up patient-years was not reported.

*Rates calculated based on data in final publication.

**In RESPECT, serious AF was adjudicated by an independent board of physicians. In REDUCE, it was determined by the local investigator.




AMPLATZER™ PFO OCCLUDER

Device Specifications

SIZING AND DEVICE SELECTION

AMPLATZER™ PFO Occluder

Model/Reorder Number	Right Atrial Disc Diameter [A]	Left Atrial Disc Diameter [B]	Minimum Recommend Sheath Size
9-PFO-018	18 mm	18 mm	8 F; 45° curve
9-PFO-025	25 mm	18 mm	8 F; 45° curve
9-PFO-030	30 mm	30 mm	8 F; 45° curve
9-PFO-035	35 mm	25 mm	9 F; 45° curve



DELIVERY SYSTEM

AMPLATZER™ TorqVue® 45° Delivery System

Model/Reorder Number	Sheath Size	Tip Angle	Sheath Inner Diameter	Sheath Outer Diameter	Usable Length
9-ITV08F45/80	8 F	45°	2.69 mm/0.11 inch	3.45 mm/0.14 inch	80 cm
9-ITV09F45/80	9 F	45°	3.00 mm/0.12 inch	3.81 mm/0.15 inch	80 cm

ANCILLARY PRODUCTS

AMPLATZER™ Guidewire

Model/Reorder Number	Diameter	Body	Tip Description	Usable Length
9-GW-002	0.035 inch	Super Stiff	1.5 mm, Modified J-tip	260 cm

REFERENCES

1. Saver JL, Carroll JD, Thaler DE, et al. Long-term outcomes of patent foramen ovale closure or medical therapy after stroke. *N Engl J Med* 2017; 377: 1022-32. 2. Tobis J, Charles A, Silbertson D, et al. Prospective, randomized investigation to evaluate incidence of headache reduction in subjects with frequent migraine and PFO using the AMPLATZER PFO occluder to medical management. *J Am Coll Cardiol* 2017; 70:2766-74. 3. Meier B, Kalesan B, Mattle HP, et al. Percutaneous closure of patent foramen ovale in cryptogenic embolism. *N Engl J Med* 2013; 368:1083 - 91. 4. Mas J-L, Derumeaux G, Guillon B, et al. Patent foramen ovale closure or anticoagulation vs. antiplatelets after stroke. *N Engl J Med* 2017;377:1011-21 and supplementary appendix. 5. Lee PH, Song J-K, Kim JS, et al. Cryptogenic Stroke and High-Risk Patent Foramen Ovale: The DEFENSE-PFO Trial, *Journal of the American College of Cardiology* (2018), doi: 10.1016/j.jacc.2018.02.046. 6. Heinrich P, Mattle, Stefan Evers, David Hildick-Smith, et al. Percutaneous closure of patent foramen ovale in migraine with aura, a randomized controlled trial, *European Heart Journal*, Volume 37, Issue 26, 7 July 2016, Pages 2029–2036. 7. Based on internal lab testing. Data on file at Abbott. 8. Søndergaard L, Kasner SE, Rhodes JF, et al. Patent foramen ovale closure or antiplatelet therapy for cryptogenic stroke. *N Engl J Med* 2017; 377:1033-42.

For more information about the AMPLATZER™ PFO Occluder or the RESPECT trial, contact your Abbott sales representative, or visit CRYPTOGENICSTROKE.COM.

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at manuals.sjm.com for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Abbott Vascular International BVBA
Park Lane, Culliganlaan 2b, 1831 Diegem, Belgium

™Indicates a trademark of the Abbott group of companies.

‡ Indicates a third party trademark, which is property of its respective owner.

©2018 Abbott. All rights reserved. SJM-AMPLP-0418-0090a | Item approved for international use only.

