


AMPLATZER Septal Occluder  
Atrial Septal Defect Closure



ST. JUDE MEDICAL  
MORE CONTROL. LESS RISK.


# AMPLATZER™ Septal Occluder

## Atrial Septal Defect Closure

The AMPLATZER Septal Occluder is a self-expanding, double-disc device designed for closure of atrial septal defects.

### Product Highlights

- Self-centering waist fills the defect
- Device can be easily recaptured and redeployed for optimal placement
- Nitinol and interwoven polyester promote occlusion and tissue in-growth
- Safe and efficacious: 97.2% closure rate at 6 months<sup>1</sup>; 1.6% major adverse events<sup>2</sup>



### Ordering Information

Contents: 1 occluder

| Model/Reorder Number | Device Size/<br>Waist Diameter<br>(mm) [A] | Waist<br>Length<br>(mm) [B] | Right Atrial<br>Disc Diameter<br>(mm) [C] | Left Atrial<br>Disc Diameter<br>(mm) [D] | Minimum<br>Recommended<br>Sheath Size* |
|----------------------|--|-----------------------------|---|--|--|
| 9-ASD-004            | 4  | 3                           | 12  | 16                                       | 6 F; 45° Curve                         |
| 9-ASD-005            | 5  | 3                           | 13  | 17                                       | 6 F; 45° Curve                         |
| 9-ASD-006            | 6  | 3                           | 14  | 18                                       | 6 F; 45° Curve                         |
| 9-ASD-007            | 7  | 3                           | 15  | 19                                       | 6 F; 45° Curve                         |
| 9-ASD-008            | 8  | 3                           | 16  | 20                                       | 6 F; 45° Curve                         |
| 9-ASD-009            | 9  | 3                           | 17  | 21                                       | 6 F; 45° Curve                         |
| 9-ASD-010            | 10   | 3                           | 18  | 22                                       | 6 F; 45° Curve                         |
| 9-ASD-011            | 11   | 4                           | 21  | 25                                       | 7 F; 45° Curve                         |
| 9-ASD-012            | 12   | 4                           | 22  | 26                                       | 7 F; 45° Curve                         |
| 9-ASD-013            | 13   | 4                           | 23  | 27                                       | 7 F; 45° Curve                         |
| 9-ASD-014            | 14   | 4                           | 24  | 28                                       | 7 F; 45° Curve                         |
| 9-ASD-015            | 15   | 4                           | 25  | 29                                       | 7 F; 45° Curve                         |
| 9-ASD-016            | 16   | 4                           | 26  | 30                                       | 7 F; 45° Curve                         |
| 9-ASD-017            | 17   | 4                           | 27  | 31                                       | 7 F; 45° Curve                         |
| 9-ASD-018            | 18   | 4                           | 28  | 32                                       | 8 F; 45° Curve                         |
| 9-ASD-019            | 19   | 4                           | 29  | 33                                       | 8 F; 45° Curve                         |
| 9-ASD-020            | 20   | 4                           | 30  | 34                                       | 9 F; 45° Curve                         |
| 9-ASD-022            | 22   | 4                           | 32  | 36                                       | 9 F; 45° Curve                         |
| 9-ASD-024            | 24   | 4                           | 34  | 38                                       | 9 F; 45° Curve                         |
| 9-ASD-026            | 26   | 4                           | 36  | 40                                       | 10 F; 45° Curve                        |
| 9-ASD-028            | 28   | 4                           | 38  | 42                                       | 10 F; 45° Curve                        |
| 9-ASD-030            | 30   | 4                           | 40  | 44                                       | 10 F; 45° Curve                        |
| 9-ASD-032            | 32   | 4                           | 42  | 46                                       | 12 F; 45° Curve                        |
| 9-ASD-034            | 34   | 4                           | 44  | 50                                       | 12 F; 45° Curve                        |
| 9-ASD-036            | 36   | 4                           | 46  | 52                                       | 12 F; 45° Curve                        |
| 9-ASD-038            | 38   | 4                           | 48  | 54                                       | 12 F; 45° Curve                        |
| 9-ASD-040**          | 40   | 4                           | 50  | 56                                       | 12 F; 45° Curve                        |

\* AMPLATZER™ TorqVue™ 45° Delivery System is recommended for use with the AMPLATZER Septal Occluder.

\*\* Device size not available for sale in the U.S.